



vW Select™

Ristocetin CoFactor Activity Assay

Catalog No. 106730

All-In-One Ristocetin CoFactor Activity Assay System

vW Select™ is an All-In-One Ristocetin CoFactor Activity Assay System comprised of reagents, control plasmas, and lot specific technical information optimized to elicit the agglutination of lyophilized platelets in the presence of the patient's Platelet Poor Plasma and Ristocetin.



- Reagents Optimized to More Effectively Perform Ristocetin CoFactor Activity Assay
- Offers a Three or Four Point Standard Curve
- Replaces Lab Developed Tests for Regulatory Compliance
- Achieves Linear Standard Curves
- Improves Test Performance, Sensitivity, and Repeatability
- Reduces Assay and Laboratory Test Variables
- Suitable for Veterinary Use

vW Select™ delivers superior repeatability in preparing standard reference curves as well as patient testing when compared to other commercial or Lab Developed Tests.

COMPLETE TEST KIT

NET CONTENTS:

Lyophilized Platelets	2 x 10.0mL
AggRecetin	2 x 15.0mg
vW Normal Reference Plasma	5 x 0.5mL
vW Abnormal Control Plasma	5 x 0.5mL
vW Normal Control Plasma	5 x 0.5mL
Tris Buffered Saline	3 x 10.0mL
AggRecetin Diluent	2 x 2.0mL



155 Gibraltar Road
Horsham, PA 19044 USA



An ISO 13485 Registered Company
www.biodatacorp.com



mdt Europa GmbH
Langenhagener Str. 71
D-30855 Langenhagen
GERMANY



Alpha Laboratories
40 Parham Drive Eastleigh
SO50 4NU Hampshire
UNITED KINGDOM



Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com