

AGG/PAK™ 5 Combo Kit

Platelet Aggregation Reagent Test Kit

Catalog No. 107650

A Comprehensive Kit with Enhanced Efficiency for Platelet Aggregation Testing

AGG/PAK™ 5 Combo Kit is a convenience kit containing a combination of routine platelet aggregation reagents used to elicit aggregation and / or agglutination responses in Platelet Rich Plasma (PRP). This Test Kit includes ADP, Arachidonic Acid, Collagen, Epinephrine, and Ristocetin Reagents.



**Conveniently Packaged for
Laboratory Efficiency**

- Evaluates Platelet Function Disorders
- Identifies von Willebrand Disease (vWD) Using Ristocetin Induced Platelet Aggregation Testing (RIPA)
- Follows ISTH / SCC Proposed Aggregation Panel
- Assesses Anti-Platelet Therapy Effectiveness
- Offers Insights into Clot Formation Dynamics
- Aids in Thrombotic Event Prevention
- Simplifies Laboratory Inventory Control

AGG/PAK™ 5 Combo Kit Reagents induce primary platelet responses, prompting shape change and aggregation. Patients with various acquired and inherited platelet function disorders show abnormal responses. A Light Transmission Platelet Aggregometer effectively captures these changes, displaying parameters such as extent and rate of aggregation. **AGG/PAK™ 5 Combo Kit** Reagents are lyophilized for long-term storage. All Reagents meet CLIA, ISTH, and other regulatory requirements, ensuring the highest standards of quality and compliance.

NET CONTENTS:

ADP	<i>Adenosine-5'-Diphosphate</i>	1 x 0.5 mL	Epinephrine	<i>Adrenaline</i>	1 x 0.5 mL
Arachidonic Acid	<i>Sodium Arachidonate</i>	1 x 0.5 mL	Ristocetin	<i>Ristocetin A Sulfate</i>	1 x 0.5 mL
Collagen	<i>Soluble Calf Skin, Type I</i>	1 x 0.5 mL			



155 Gibraltar Road
Horsham, PA 19044 USA

REF
107650



An ISO 13485 Registered Company
www.biodatacorp.com



mdi Europa GmbH
Langenhagener Str. 71
D-30855 Langenhagen
GERMANY



Alpha Laboratories
40 Parham Drive Eastleigh
S050 4NU Hampshire
UNITED KINGDOM



Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com